

PRESENTATION OF THE BLIFT + PROGRAMS

Treatment family	Pathologies	Accessories	Modes	stages	Tps/ Pressure
SCARS	Retractable scars	Suction cup	Pulsed	Mobilize	5 min/ 500mbar
		Suction cup	Continuous	Smooth	5 min/ 530mbar
	Fibrotic scars	Suction cup	Pulsed	Mobilize	5 min/ 400mbar
		Suction cup	Continuous	Smooth	5 min/ 400mbar
	Adherent scars	Suction cup	Pulsed	Mobilize	5 min/ 500mbar
		Suction cup	Continuous	Smooth	5 min/ 600mbar
	Cesarean section	R35 suction cup	Continuous	Vasculariser	5 min/ 200mbar
		Palpate and roll	Continuous	Smooth	10min/ 300mbar
DRAINAGE	Edema	Suction cup	Pulsed	Mobilize	5 min/ 530mbar
		Suction cup	Continuous	Drain	10min/ 530mbar
	Upper limbs	Suction cup	Pulsed	Stimulate	5 min/ 530mbar
		Suction cup	Pulsed	Mobilize	5 min/ 530mbar
	Lower limbs	Suction cup	Continuous	Drain	10min/ 530mbar
		Suction cup	Pulsed	Stimulate	5 min/ 530mbar
		Suction cup	Pulsed	Mobilize	5 min/ 530mbar
		Suction cup	Continuous	Drain	10min/ 530mbar
BURNS	Suction cup / Palpate and roll	Pulsed	Mobilize	5 min/ 150mbar	
	Suction cup	Continuous	Soften	5 min/ 200mbar	
TENDINOPATHY	Suction cup	Pulsed	Relieve	5 min/ 300mbar	
	Suction cup	Continuous	Mobilize	5 min/ 300mbar	

PRESENTATION OF THE BLIFT + PROGRAMS

Treatment family	Pathologies	Accessories	Modes	stages	Tps/ Pressure
RACHIS	Cervicalgie	R15 / R25 suction cup	Continuous	Vasculariser	5 min/ 200mbar
		Suction cup	Pulsed	Mobilize	5 min/ 200mbar
		Suction cup	Continuous	Soften	5 min/ 200mbar
	Dorsalgia	Suction cup	Continuous	Vasculariser	5 min/ 500mbar
		Suction cup	Pulsed	Mobilize	5 min/ 500mbar
		Suction cup	Continuous	Soften	5 min/ 500mbar
	Low back pain	Suction cup	Continuous	Vasculariser	5 min/ 200mbar
		Suction cup	Pulsed	Mobilize	5 min/ 200mbar
		Suction cup	Continuous	Soften	5 min/ 200mbar
	Sciatalgie	Suction cup	Continuous	Vasculariser	5 min/ 200mbar
		Suction cup	Pulsed	Mobilize	5 min/ 200mbar
		Suction cup	Continuous	Soften	5 min/ 200mbar
Irradiated skins	Suction cup	Continuous	Mobilize	5 min/ 200mbar	
	Suction cup	Continuous	Soften	5 min/ 300mbar	



Mazet Electronics ZA Route de Tence, 43520 le Mazet-Saint-Voy